

# BIO-RAD Lyphocheck Immunoassay Plus Control

For the control of cTSH, cPROG and cCOR

**FOR VETERINARY USE ONLY!**



**Item code 105382**

**LOT 40412**

**Expiry date 30.09.2024**

## Handling

To be reconstituted with 5 ml of distilled water, mix the control gently and wait for 30min.

## Target values conventional units

Element RC3X						
Analyte	Unit	Mean	Range	Unit	Mean	Range
cTSH	ng/ml	1,0	0,8 – 1,2	ng/ml	1,0	0,8 – 1,2
cPROG	ng/ml	10,5	8,4 – 12,6	nmol/L	33,4	26,7 – 40,1
cCOR	µg/dl	5,6	4,4 – 6,8	nmol/L	154,6	121,4 – 187,7

\*\* QC mode is using conventional units only

## Artikel – Nr. 105382

### Fuji DRI-CHEM Immuno AU 10V Qualitätskontrolle TSH, Cortisol, Progesteron

Lyphocheck Immunoassay Plus Control Level 2

HANDLING: Mit 5 ml destilliertem Wasser auflösen und für **30** Minuten bei Raumtemperatur inkubieren

Lot **40412**  
Verfall **30.09.2024**

Parameter	Zielwert	Konv. Einheit Zielbereich	Einheit	Zielwert	SI-Einheit Zielbereich	Einheit
TSH	2,15	1,72 – 2,58	ng/mL	2,15	1,72 – 2,58	ng/mL
COR (DIL) LOT 31042 or earlier	9,6	7,7 – 11,5	µg/dL	264,7	212,52 – 317,4	nmol/L
COR (DIL) LOT 31043 or later	14,8	11,8 – 17,8	µg/dL	408,5	325,7 – 491,3	nmol/L
PRG	23,69	18,95 – 28,43	ng/mL	76,32	60,42 – 89,0	nmol/L

