

# Ergebnisprotokoll FDC IMMUNO AU10V

Serum-Kontrollwerte  
Kontrolle für T4, TSH, COR, PRG, BA  
**FOR VETERINARY USE ONLY!**



## Artikel – Nr 105379 Fuji DRI-CHEM Immuno AU 10V Qualitätskontrolle T4

Lyphocheck Immunoassay Plus Control Level 1  
HANDLING: Mit 5 ml destilliertem Wasser auflösen und für 15 Minuten bei Raumtemperatur inkubieren

Lot 40421  
Verfall 31.03.2025

Parameter	Konv. Einheit			SI-Einheit		
	Zielwert	Zielbereich	Einheit	Zielwert	Zielbereich	Einheit
T4	3,59	2,87 – 4,31	µg/dL	46,2	36,9 – 55,5	nmol/L

## Artikel – Nr 105382 Fuji DRI-CHEM Immuno AU 10V Qualitätskontrolle TSH, Cortisol, Progesteron

Lyphocheck Immunoassay Plus Control Level 2  
HANDLING: Mit 5 ml destilliertem Wasser auflösen und für 15 Minuten bei Raumtemperatur inkubieren

Lot 40422  
Verfall 31.03.2025

Parameter	Konv. Einheit			SI-Einheit		
	Zielwert	Zielbereich	Einheit	Zielwert	Zielbereich	Einheit
TSH	2,29	1,83 – 2,74	ng/mL	2,29	1,83 – 2,74	ng/mL
COR (DIL) LOT 31042 or earlier	9,0	7,2 – 10,8	µg/dL	248,3	198,6 – 297,9	nmol/L
COR (DIL) LOT 31043 or later	14,9	11,9 – 17,9	µg/dL	411,0	328,3 – 493,8	nmol/L
PRG	23,71	18,97 – 28,45	ng/mL	75,40	60,32 – 90,47	nmol/L

## Artikel – Nr 105383 Fuji DRI-CHEM Immuno AU 10V Qualitätskontrolle Gallensäuren Level 2

Total Bile Acids Assay Control Level 2  
HANDLING: Mit 1 ml destilliertem Wasser auflösen und für 30 Minuten bei Raumtemperatur inkubieren

Lot BAC1200-01-06  
Verfall 04.05.2022

Parameter	Konv. Einheit			SI Einheit		
	Zielwert	Zielbereich	Einheit	Zielwert	Zielbereich	Einheit
BA	23,3	18,7 – 28,0	µmol/L	23,3	18,7 – 28,0	µmol/L

Wenn die Ergebnisse außerhalb des Zielbereichs liegen, wenden Sie sich bitte an Ihren Technischen Service.



# QC value sheet FDC IMMUNO AU10V

Serum Control values  
For the control of T4, TSH, COR, PRG, BA  
**FOR VETERINARY USE ONLY!**



## Article – No 105379 Fuji DRI-CHEM Immuno AU 10V Quality control T4

Lyphocheck Immunoassay Plus Control Level 1  
HANDLING: To be reconstituted with 5 ml of distilled water for 15min at room temperature

Lot 40411  
Exp. 30.09.2022

Analyte	Conv. Unit			International Unit		
	Mean	Range	Units	Mean	Range	Units
T4	3,59	2,87 – 4,31	µg/dL	46,2	36,9 – 55,5	nmol/L

## Article – No 105382 Fuji DRI-CHEM Immuno AU 10V Quality control TSH, Cortisol, Progesteron

Lyphocheck Immunoassay Plus Control Level 2  
HANDLING: To be reconstituted with 5 ml of distilled water for 15min at room temperature

Lot 40411  
Exp. 30.09.2022

Analyte	Conv. Unit			International Unit		
	Mean	Range	Units	Mean	Range	Units
TSH	2,29	1,83 – 2,74	ng/mL	2,29	1,83 – 2,74	ng/mL
COR (DIL) LOT 31042 or earlier	9,0	7,2 – 10,8	µg/dL	248,3	198,6 – 297,9	nmol/L
COR (DIL) LOT 31043 or later	14,9	11,9 – 17,9	µg/dL	411,0	328,3 – 493,8	nmol/L
PRG	23,71	18,97 – 28,45	ng/mL	75,40	60,32 – 90,47	nmol/L

## Article – No 105383 Fuji DRI-CHEM Immuno AU 10V Quality Bile Acids Level 2

Total Bile Acids Assay Control Level 2  
HANDLING: To be reconstituted with 1 ml of distilled water for 30min at room temperature

Lot BAC1200-01-06  
Exp. 04.05.2022

Analyte	Conv. Unit			International Unit		
	Mean	Range	Units	Mean	Range	Units
BA	23,3	18,7 – 28,0	µmol/L	23,3	18,7 – 28,0	µmol/L

If the results are out of range, please contact your Technical Service.

