


# Validation of a new canine species-specific C-reactive protein assay on the Pentra 400

S.Klenner<sup>1</sup>, S. Zielinsky<sup>2</sup>, N. Kneier<sup>1</sup>, J. Klotz<sup>1</sup>, N. Bauer<sup>2</sup>, A. Moritz<sup>2</sup>

<sup>1</sup>scil animal care company GmbH, Viernheim, <sup>2</sup>Department of Veterinary Clinical Sciences, Clinical Pathology and Clinical Pathophysiology, Justus-Liebig-University Giessen, Germany



## Introduction & Aim

Rapid immunoturbidimetric detection of canine C-reactive protein (CRP) replaced CRP detection with ELISA-methodology but was hampered by the use of a human CRP test and labour intensive by preparation of pool plasma as canine controls. Costly purified canine CRP was further needed to establish a species-specific calibration curve.

Now a new canine-specific immunoturbidimetric CRP assay is available which provides canine-specific reagents, controls and calibrators.

Aim of the study was the validation of this new assay (by **Gentian AS**, Moss, Norway) against the previously validated immunoturbidimetric assay from **Randox Laboratories Ltd** (Crumlin, UK) [1] on the wet chemistry analyzer **Pentra 400** (ABX horiba diagnostics).

## Material and Methods

### Inaccuracy: Linearity under dilution

- Patient serum 281.3 mg/l CRP diluted with NaCl to 7 different levels, each level was measured in triplicates

### Intra- and inter-assay precision:

- 4 different levels of patient sera:  
A: 7.15 mg/l, B: 58.36 mg/l, C: 103.86 mg/l, D: 272.07 mg/l CRP
- measured daily, twice in the morning and twice in the evening for 5 consecutive days

### Interference:

- aliquots of canine patient serum of 35,5 mg/L CRP spiked with 800 mg/l bilirubine, 5 g/l hemoglobin, or 10 g/l intralipid were analyzed in triplicates and compared to aliquots spiked with either NaOH, NaCl or A. dest.

### Evaluation of possible prozone effect:

- Three different CRP samples were used to assess hook-effect: 455 mg/l, 676 mg/l, 890 mg/l

### Method comparison:

- 279 sera of healthy and diseased dogs
- Comparison with immunoturbidimetric assay from Randox Laboratories Ltd, Crumlin, UK [1]

## Statistical analysis

- Mean and standard deviation (SD), coefficient of variation (CV)
- Spearman's rank correlation coefficient
- Bland-Altman analysis
- Quality specifications were set according to the ASVCP TE<sub>a</sub> guidelines [2]:
  - TE<sub>a</sub> optimal/ desired: 14.79%/ 29.58%
  - CV optimal/ desired: 6.06%/ 12.16%
  - bias optimal/ desired: 4.76%/ 9.52 %

## Conclusions

The new Gentian canine CRP assay precisely and accurately detects canine CRP in comparison to the previously used immunoturbidimetric method. Advantage is the availability of a whole test kit including species-specific reagents, calibrators, and controls.

## References

- [1] Klenner S, Bauer N, Moritz A. Evaluation of three automated human immunoturbidimetric assays for the detection of C-reactive protein in dogs. J Vet Diagn Invest. 2010 Jul;22(4):544-52.  
[2] K.E. Harr, B. Flatland, M. Nabity, K.Freeman, ASVCP Allowable Total Error Recommendations for Biochemistry, <http://www.asvcp.org/pubs/qas/index.cfm>, assessed 24<sup>th</sup> of September, 2013

The Eurostars Programme is powered by EUREKA and the European Community


This study was performed as part of the project „MAMMALIAN CRP“ within the framework of the Eurostars Programme and supported by the German Federal Ministry of Education and Research (FKZ 01QE1110).  
The sole responsibility for the study's content lies with the authors.



## Results

### Inaccuracy:

Linearity under dilution was present up to 280 mg/l.



### Precision:

Optimal quality specifications could be met. CV of the lowest CRP level was lower than the desired CV% outlined by the ASVCP.

Level	Within-run [mg/l]			Between-run [mg/l]		
	Mean	SD	CV%	Mean	SD	CV%
A	7.15	0.87	12.12	7.15	0.56	7.84
B	58.36	1.08	1.85	58.36	1.25	2.14
C	103.86	0.88	0.85	103.86	0.91	0.88
D	272.07	1.85	0.68	272.07	4.98	1.83

### Interference:

Mean absolute bias calculated was 0.57 mg/l, 0.13 mg/l, and 1.63 mg/l for the interferent bilirubine, hemoglobin and intralipid. At a clinical decision level of i.e. 30 mg/l CRP, the optimal TE<sub>a</sub> would be 4.44mg/l in absolute numbers. Observed systematic error is lower than the absolute total allowable error in all interfering substances.

### Evaluation of possible prozone effect:

Recovery rates for the detection of 455 mg/l CRP were 130% for 676 mg/l 121% and for 890 mg/l 28%. No prozone effect was present up to a concentration of 676 mg/l CRP.

### Method Comparison:

- Spearman's rank correlation coefficient ( $r$ ),  $n=263$ , 12 excluded due to missing data,  $r= .99$
- Bland-Altman analysis,  $n=263$ , 12 excluded due to missing data, bias -1.2%

