



Parameter	Comprehensive Plus	Pre-Surgical	Electrolytes	Liver	Equine Plus / Large Animal	Kidney	Bovine Monitoring	Kidney Monitoring	Avian Reptilian	General Health	Glucose
ALB	X	X		X	X	X			X	X	
ALP	X	X		X	X					X	
ALT	X	X		X						X	
AMY	X					X					
AST				X	X		X		X	X	
BUN	X	X			X	X		X		X	
CA	X		X		X	X	X		X	X	
TC				X						X	
CK					X		X		X		
CL			X								
TCO ₂			X			X					
CREA	X	X			X	X		X		X	
GGT				X	X					X	
GLU	X	X			X	X			X	X	X
K	X		X		X				X	X	
LDH					X						
NA	X		X		X				X	X	
MG			X		X		X				
PHOS	X		X		X	X	X	X	X	X	
TB	X			X	X					X	
TBA				X					X		
TG					X						
TP	X	X		X	X				X	X	
URIC ACID									X		
BUN/CREA*	X	X			X	X		X		X	
GLOB*	X	X		X	X				X	X	
A/G*	X	X		X	X				X	X	

*berechnete Parameter